

Invia® Liberty™ and Foam Dressing Kits



SIMPLIFYING NEGATIVE PRESSURE WOUND THERAPY



For User-Friendly Negative Pressure Wound Therapy


Invia® Liberty™ and Avance® Foam Dressing Kits
Providing a flexible solution for Negative Pressure Wound Therapy

Negative Pressure Wound Therapy (NPWT) has become a preferred choice of care for treating acute and chronic wounds. Invia® Liberty™ delivers NPWT that satisfies the demands of clinicians and patients.

Invia Liberty's innovative design delivers excellent functionality with its simple set-up and easy operation. Invia Liberty can be used both at hospital and at home, easing transition and avoiding extra instruction time for patients.

Patients appreciate the ease-of-use, quiet operation and compact size of the Liberty. The device can be carried in a discreet bag, facilitating mobility. The docking station allows for easy recharging of the Invia Liberty.

Therapy data can be checked on the built-in LCD screen or stored and uploaded to a PC.



- (a) Disposable canister 300 ml & 800 ml options with solidifier
- (b) One-touch canister release button
- (c) Easy to read LCD screen for intuitive operation and graphic capability
- (d) Double lumen tubing (PVC) with integrated protection filter and tube clamp
- (e) Covered, waterspray safe power adapter port
- (f) Covered, waterspray safe USB port for software updating and therapy history upload to a PC

Avance® Foam Dressing Kits

- Promotes formation of granulation tissues in wounds
- Multiple sizes provide flexibility to manage different types of wounds
- Bleeding and exudate can be more easily monitored when using green foam¹
- Hydrophobic material assists with removal of exudate
- High tensile strength with minimal fragmenting when cutting to wound size²
- No difference has been found in the quantity or characteristics of the granulation tissue formed under the green or black foam¹



¹Malmjö et al. *J of Wound Care* 2011 Jun;20(6):294-9. Green foam, black foam or gauze for NPWT: effects on granulation tissue formation.

²Data on file with Molnlycke Healthcare.

Small Dressing Kit

- 1 Avance® Foam (8 x 10 x 3 cm)
- 1 Avance Transparent Film (20 x 27cm)
- 1 Avance Transfer Pad (ESI)

Medium Dressing Kit

- 1 Avance Foam (12 x 20 x 3cm)
- 2 Avance Transparent Film (20 x 27cm)
- 1 Avance Transfer Pad (ESI)

Large Dressing Kit

- 1 Avance Foam (15 x 26 x 3cm)
- 2 Avance Transparent Film (20 x 27cm)
- 1 Avance Transfer Pad (ESI)

Medical Vacuum Technology
 for Healthcare Professionals

Please contact us or your local Medela representative for details.

Federal law restricts this device to sale or rental by or on the order of a physician. Avance Dressing Kits distributed by Medela, Inc.

Medela and Invia are registered trademarks of Medela Holding AG. Liberty is a trademark of Medela, Inc. Avance is a registered trademark of Molnlycke.

Medela AG
 Lättichstrasse 4b
 6341 Baar / Switzerland
 www.medela.com

Medela, Inc.
 1101 Corporate Drive
 McHenry, IL 60050
 USA
 Phone +1 877 735 1626
 Fax +1 815 307 8942
 suction@medela.com
 www.medelasuction.com

Local contact:



Expertise You Trust • Service You Deserve
 Phone 888.573.1400 | Fax 888.665.4199